

Anti-CD85j/ILT2 antibody {ABflo® 488} (STJ11108536)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | ABflo® 488 |
| Purification | Affinity purification |
| Dilution Range | FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume |
| Formulation | PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype | IgG |
| Molecular Weight | 49kDa/70kDa/71kDa |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 10859 |
| Gene Symbol | LILRB1 |
| UniProt ID | LIRB1_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|