

Anti-MMP9 antibody {ABflo® 647} (STJ11108546)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC (intra) , 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	78kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	4318
Gene Symbol	MMP9
UniProt ID	MMP9_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---