

Anti-CD305/LAIR1 antibody {ABflo® 647} (STJ11108585)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	23kDa/29kDa/31kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.
TARGET INFORMATION	
Gene ID	3903
Gene Symbol	LAIR1
UniProt ID	LAIR1_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.