

Anti-GP6/GPVI antibody {ABflo® 647} (STJ11108593)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	FC
	Host / Source	Rabbit
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	ABflo® 647
	Purification	Affinity purification
	Dilution Range	FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume
	Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
	Isotype	IgG
	Molecular Weight	37kDa
	Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.
TARGET INFORMATION		
	Gene ID	51206
	Gene Symbol	GP6
	UniProt ID	GPVI_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.