

Anti-Cytokeratin 19 (KRT19) antibody {ABflo® 594} (STJ11108630)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	ABflo® 594
Purification	Affinity purification
Dilution Range	FC (intra) , 5 Mu l per 10^6 cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	44kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.
TARGET INFORMATION	
Gene ID	3880
Gene Symbol	KRT19
UniProt ID	K1C19_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.