

Anti-IL-17RB antibody {ABflo® 488} (STJ11108651)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	31kDa/55kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	55540
Gene Symbol	IL17RB
UniProt ID	I17RB_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---