

**Anti-CD177 antibody {ABflo® 647} (STJ11108674)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                 |
| Conjugation         | ABflo® 647                                                                                                                 |
| Purification        | Affinity purification                                                                                                      |
| Dilution Range      | FC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume                                                                    |
| Formulation         | PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype             | IgG                                                                                                                        |
| Molecular Weight    | 15kDa/27kDa/46kDa                                                                                                          |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze.                                                     |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">57126</a>       |
| Gene Symbol | <a href="#">CD177</a>       |
| UniProt ID  | <a href="#">CD177_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|