

Anti-CD127/IL-7R Alpha antibody {ABflo® 647} (STJ11108725)**GENERAL INFORMATION**

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rabbit |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | ABflo® 647 |
| Purification | Affinity purification |
| Dilution Range | FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume |
| Formulation | PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype | IgG |
| Molecular Weight | 51kDa |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 16197 |
| Gene Symbol | IL7r |
| UniProt ID | IL7RA_MOUSE |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|