

**Anti-BTN3A1/2/3 antibody {ABflo® 488} (STJ11108733)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	58kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

**TARGET INFORMATION**

Gene ID [10384](#)  
[11118](#)  
[11119](#)  
[BTN3A3](#)  
[BTN3A2](#)  
[BT3A3\\_HUMAN](#)  
[BT3A2\\_HUMAN](#)  
[BT3A1\\_HUMAN](#)

**ADDITIONAL INFORMATION**

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.