

**Anti-CD101/IGSF2 antibody {ABflo® 488} (STJ11108748)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IF/ICC/FC
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	IF/ICC, 1:50-1:200IFC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	115kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">9398</a>
Gene Symbol	<a href="#">CD101</a>
UniProt ID	<a href="#">IGSF2_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---