

## Anti-CD8b antibody {ABflo® 647} (STJ11108764)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | ABflo® 647   |
| Purification        | Affinity purification  |
| Dilution Range      | FC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume  |
| Formulation         | PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype             | IgG  |
| Molecular Weight    | 24kDa  |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze.   |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">926</a>        |
| Gene Symbol | <a href="#">CD8B</a>       |
| UniProt ID  | <a href="#">CD8B_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|