

Anti-CD235a/Glycophorin A antibody {ABflo® 594} (STJ11108782)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Conjugation | ABflo® 594 |
| Purification | Affinity purification |
| Dilution Range | FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume |
| Formulation | PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3. |
| Isotype | IgG |
| Molecular Weight | 13kDa/16kDa |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 2993 |
| Gene Symbol | GYPA |
| UniProt ID | GLPA_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|