

**Anti-CD126/IL-6R Alpha chain antibody {ABflo® 594} (STJ11108791)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	ABflo® 594
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	50kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

**TARGET INFORMATION**

Gene ID	<a href="#">16194</a>
Gene Symbol	<a href="#">IL6ra</a>
UniProt ID	<a href="#">IL6RA_MOUSE</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---