

## Anti-CD24 antibody {ABflo® 594} (STJ11108823)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/FC
Host / Source	Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 594
Purification	Affinity purification
Dilution Range	IF/ICC, 1:50-1:200IFC, 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	8kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">12484</a>
Gene Symbol	<a href="#">Cd24</a>
UniProt ID	<a href="#">CD24_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---