

Anti-CD6 antibody {ABflo® 647} (STJ11108838)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	64kDa/68kDa/72kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.
TARGET INFORMATION	
Gene ID	12511
Gene Symbol	Cd6
UniProt ID	CD6_MOUSE
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.