

Anti-Rabbit IgG antibody {ABflo® 700} (STJ11108892)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Rabbit

PRODUCT PROPERTIES

Conjugation	ABflo® 700
Purification	Affinity purification
Dilution Range	FC, 5 µl per 10 ⁶ cells in 100 µl volume
Formulation	PBS with 0.03% proclin300, 0.2% BSA, pH7.3.
Isotype	IgG
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---