

## Anti-IL-9 antibody {ABflo® 647} (STJ11108929)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC (intra) , 5 Mu l per 10 <sup>6</sup> cells in 100 Mu l volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	16kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">16198</a>
Gene Symbol	<a href="#">IL9</a>
UniProt ID	<a href="#">IL9_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------