

Anti-CD94 antibody {ABflo® 647} (STJ11108959)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC, 5 µl per 10 ⁶ cells in 100 µl volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	21kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	16643
Gene Symbol	Klrd1
UniProt ID	KLRD1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---