

Anti-CD337/NKp30 antibody {ABflo® 450} (STJ11108982)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 450
Purification	Affinity purification
Dilution Range	FC, 5 µl per 10 ⁶ cells in 100 µl volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	16kDa/17kDa/18kDa/19kDa/20kDa/21kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	259197
Gene Symbol	NCR3
UniProt ID	NCTR3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---