

Anti-CD87/PLAUR antibody {PE} (STJ11108996)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Purification	Affinity purification
Dilution Range	FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	37kDa/31kDa/32kDa/36kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	5329
Gene Symbol	PLAUR
UniProt ID	UPAR_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---