

Anti-Integrin alpha 1/CD49a antibody {ABflo® 488} (STJ11109001)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 488
Purification	Affinity purification
Dilution Range	FC, 5 µl per 10 ⁶ cells in 100 µl volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	130kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	109700
Gene Symbol	Itga1
UniProt ID	ITA1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---