

Anti-CD94 antibody {PE} (STJ11109010)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | PE |
| Purification | Affinity purification |
| Dilution Range | FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume |
| Formulation | PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3. |
| Isotype | IgG |
| Molecular Weight | 17kDa/20kDa |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 3824 |
| Gene Symbol | KLRD1 |
| UniProt ID | KLRD1_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|