

Anti-CD200/OX2 antibody {PE} (STJ11109055)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rabbit |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | PE |
| Purification | Affinity purification |
| Dilution Range | FC, 5 Mu l per 10 ⁶ cells in 100 Mu l volume |
| Formulation | PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3. |
| Isotype | IgG |
| Molecular Weight | 31kDa |
| Storage Instruction | Store at 4°C for up to 1 year from the date of receipt. Do not freeze. |

TARGET INFORMATION

| | |
|-------------|--|
| Gene Symbol | Cd200 Uniprot ID= OX2G_MOUSE " |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|