

## Anti-CD31 antibody {Biotin} (STJ11109108)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Biotin
Purification	Affinity purification
Dilution Range	FC, $\leq 0.25 \mu\text{g}$ per million cells in 100 $\mu\text{l}$ volume
Formulation	PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.
Isotype	IgG
Molecular Weight	69kDa/77kDa/80kDa/81kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

### TARGET INFORMATION

Gene ID	<a href="#">18613</a>
Gene Symbol	<a href="#">Pecam1</a>
UniProt ID	<a href="#">PECA1_MOUSE</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---