

Anti-Phospho-PI3KP85 Alpha-Y467+P85 Beta-Y464+P55 Gamma-Y199 antibody {ABflo® 647} (STJ11109142)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	ABflo® 647
Purification	Affinity purification
Dilution Range	FC (intra) , 5 Mu l per 10 ⁶ cells in 100 Mu l volume
Formulation	PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis) , pH 7.3.
Isotype	IgG
Molecular Weight	42-54kDa/81-84kDa
Storage Instruction	Store at 4°C for up to 1 year from the date of receipt. Do not freeze.

TARGET INFORMATION

Gene ID	8503 5295 5296 PIK3R3 PIK3R1 P55G_HUMAN P85A_HUMAN P85B_HUMAN
---------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---