

## Anti-SATB2 antibody [SRM14416] (STJ11109146)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF-P/IHC-P/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:10000 IF-P, 1:200-1:800 IHC-P, 1:20000-1:40000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on
Isotype	IgG
Molecular Weight	Protein Weight 15kDa Observed Mw 15kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">23314</a>
Gene Symbol	<a href="#">SATB2</a>
UniProt ID	<a href="#">SATB2_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---