

Anti-PD-L1/CD274 antibody [SRM14418] (STJ11109148)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/ELISA/Inhibition
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	≥95% purity by SDS-PAGE
Dilution Range	FC, 1:50-1:200 Elisa, 1:500-1:2000
Isotype	IgG
Molecular Weight	Protein Weight 54kDa Observed Mw 70kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---