

Anti-EIF4G2/p97 antibody [SRM14441] (STJ11109171)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC-P/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:20000 IHC-P, 1:500-1:1000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-600 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requi
Isotype	IgG
Molecular Weight	Protein Weight 32kDa Observed Mw 35kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1982
Gene Symbol	EIF4G2
UniProt ID	IF4G2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---