

Anti-DiMethyl-Histone H3-K36 antibody [SRM14484] (STJ11109214)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/DB/IHC-P/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Other
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:10000 DB, 1:2000-1:10000 IHC-P, 1:1000-1:5000 IF/ICC, 1:50-1:200 IP, 2 Mu g-6 Mu g antibody for 300 Mu g-600 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the con
Isotype	IgG
Molecular Weight	Protein Weight 72kDa Observed Mw 130kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	8290 8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
Gene Symbol	H3-4 H3C1.H3C2.H3C3.H3C4.H3C6.
UniProt ID	H31T_HUMAN H31_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.