

## Anti-MonoMethyl-Histone H3-K27 antibody [SRM14499] (STJ11109229)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/DB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/A.thaliana/Rice/Other

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:30000 DB, 1:10000-1:30000 IF/CC, 1:50-1:200 ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 43kDa Observed Mw 43kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">8350/8351/8352/8353/8354/8355/8356/8357/8358/89688290</a>
Gene Symbol	<a href="#">H3C1.H3C2.H3C3.H3C4.H3C6.H3C7.H3C8.H3C10.H3C11.H3C12</a> <a href="#">H31_HUMAN</a> <a href="#">H31T_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---