

Anti-Pan-Succinyl Lysine antibody [SRM14508] (STJ11109238)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/DB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:5000 DB, 1:500-1:2000 IP, 0.5 μ g-4 μ g antibody for 200 μ g-600 μ g extracts of whole cells ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requireme
Isotype	IgG
Molecular Weight	Protein Weight 17kDa/20kDa Observed Mw 21kDa/24kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---