

Anti-BACH1 antibody [SRM14510] (STJ11109240)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA/ChIP
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:8000 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements. ChIP, 5 μ g antibody for 10 μ g-15 μ g of Chromatin
Isotype	IgG
Molecular Weight	Protein Weight 129kDa Observed Mw 180kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	571
Gene Symbol	BACH1
UniProt ID	BACH1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---