

Anti-DNMT1 antibody [SRM14534] (STJ11109264)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:1000-1:4000 IHC-P, 1:1000-1:4000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requi
	Isotype	IgG
	Molecular Weight	Protein Weight 48kDa Observed Mw 48kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	1786
	Gene Symbol	DNMT1
	UniProt ID	DNMT1_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.