

## Anti-ELN antibody [SRM14567] (STJ11109297)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB/IF-P/IHC-P/ELISA |
| Host / Source | Rabbit              |
| Reactivity    | Human/Rat           |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Affinity purification   |
| Dilution Range      | WB, 1:500-1:1000<br>IF-P, 1:100-1:400<br>IHC-P, 1:500-1:2000<br>ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements. |
| Isotype             | IgG   |
| Molecular Weight    | Protein Weight 102 kDa<br>Observed Mw 160 kDa   |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.  |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">2006</a>      |
| Gene Symbol | <a href="#">ELN</a>       |
| UniProt ID  | <a href="#">ELN_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|