

Anti-Galectin 3/LGALS3 antibody [SRM14575] (STJ11109305)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IF/ICC/IP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:10000-1:40000 IF/ICC, 1:200-1:800 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requ
	Isotype	IgG
	Molecular Weight	Protein Weight 22kDa Observed Mw 27 kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3958
	Gene Symbol	LGALS3
	UniProt ID	LEG3_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.