

## Anti-ALOX15 antibody [SRM14578] (STJ11109308)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A purified
Dilution Range	WB, 1:1000-1:3000 IHC-P, 1:4000-1:16000 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Isotype	IgG
Molecular Weight	Protein Weight 12kDa Observed Mw 12kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">246</a>
Gene Symbol	<a href="#">ALOX15</a>
UniProt ID	<a href="#">LOX15_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---