

## Anti-TFF1/pS2 antibody [SRM14581] (STJ11109311)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 IP, 0.5 $\mu$ g-4 $\mu$ g antibody for 200 $\mu$ g-400 $\mu$ g extracts of whole cells IF-P, 1:50-1:200 ELISA, Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements
Isotype	IgG
Molecular Weight	Protein Weight 65kDa/61kDa/66kDa/47kDa Observed Mw 68kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7031</a>
Gene Symbol	<a href="#">TFF1</a>
UniProt ID	<a href="#">TFF1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---