

Anti-VDR antibody [SRM14621] (STJ11109351)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF/ICC/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:1000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF/ICC, 1:100-1:400 IF-P, 1:100-1:400 IHC-P, 1:500-1:2000 ChIP, 5 Mu g antibody for 5 Mu g-10 Mu g of Chromatin ELISA, Recommended sta
Isotype	IgG
Molecular Weight	Protein Weight 82kDa Observed Mw 83kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7421
Gene Symbol	VDR
UniProt ID	VDR_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---