

Anti-SMARCA2/BRM antibody [SRM14623] (STJ11109353)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/IP/ChIP/ELISA
Host / Source	Rabbit
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:500-1:2000 IHC-P, 1:200-1:800 IF/ICC, 1:50-1:200 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ChIP, 5 Mu g antibody for 10 Mu g-20 Mu g of Chromatin ELISA, Recommended starting concentration is
Isotype	IgG
Molecular Weight	Protein Weight 70kDa Observed Mw 70kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	6595
Gene Symbol	SMARCA2
UniProt ID	SMCA2_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.