

Anti-CD162/PSGL-1 antibody [SRM14627] (STJ11109357)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IF/ICC/ELISA
	Host / Source	Rabbit
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:100000-1:400000 IHC-P, 1:1000-1:5000 IF/ICC, 1:100-1:500 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 61kDa Observed Mw 61kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	6404
	Gene Symbol	SELPLG
	UniProt ID	SELPL_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.