

Anti-USF1 antibody [SRM14644] (STJ11109374)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:10000 IF-P, 1:50-1:200 IHC-P, 1:500-1:1000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells
Isotype	IgG
Molecular Weight	Protein Weight 57kDa Observed Mw 57kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7391
Gene Symbol	USF1
UniProt ID	USF1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---