

Anti-CD45 antibody [SRM14649] (STJ11109379)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/mIHC/ELISA
	Host / Source	Rabbit
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB, 1:100-1:500 IHC-P, 1:5000-1:20000 mIHC, 1:5000-1:20000 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Isotype	IgG
	Molecular Weight	Protein Weight 24kDa Observed Mw 28kDa
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	5788
	Gene Symbol	PTPRC
	UniProt ID	PTPRC_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.