

Anti-eIF4A1 antibody [SRM14653] (STJ11109383)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:5000-1:20000 IHC-P, 1:400-1:2500 IF/ICC, 1:1000-1:4000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on
Isotype	IgG
Molecular Weight	Protein Weight 75kDa Observed Mw 75kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1973
Gene Symbol	EIF4A1
UniProt ID	IF4A1_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.