

## Anti-BCAT2 antibody [SRM14683] (STJ11109413)

### GENERAL INFORMATION

|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB/IF-P/IHC-P/ELISA |
| Host / Source | Rabbit              |
| Reactivity    | Human/Mouse/Rat     |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                       |
| Conjugation         | Unconjugated                                                                                                                                                                                     |
| Purification        | Affinity purification                                                                                                                                                                            |
| Dilution Range      | WB, 1:1000-1:5000<br>IF-P, 1:50-1:200<br>IHC-P, 1:50-1:200<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Isotype             | IgG                                                                                                                                                                                              |
| Molecular Weight    | Protein Weight 37kDa<br>Observed Mw 37kDa                                                                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                   |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">587</a>         |
| Gene Symbol | <a href="#">BCAT2</a>       |
| UniProt ID  | <a href="#">BCAT2_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|