

Anti-Beta III-Tubulin antibody [SRM14712] (STJ11109442)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:3000-1:20000 IP, 0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells IF-P, 1:50-1:200 IHC-P, 1:50-1:200 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your s
Isotype	IgG
Molecular Weight	Protein Weight 308kDa Observed Mw Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10381
Gene Symbol	TUBB3
UniProt ID	TBB3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---