

Anti-Histone H2B (testis specific) antibody [SRM14729] (STJ11109459)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IP/IF-P/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Other
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:6000 IP, 0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells IF-P, 1:200-1:800 IHC-P, 1:200-1:2000 ChIP, 5 Mu g antibody for 5 Mu g-15 Mu g of Chromatin ELISA, Recommended starting concentration is
Isotype	IgG
Molecular Weight	Protein Weight 11kDa Observed Mw
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	255626
Gene Symbol	H2BC1
UniProt ID	H2B1A_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.