

Anti-CD101/IGSF2 antibody [SRM14750] (STJ11109480)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IF/ICC/FC
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:6000 IF/ICC, 1:200-1:800 FC, 1:100-1:500
Isotype	IgG
Molecular Weight	Protein Weight 36kDa Observed Mw 46 kDa (LAP) /48 kDa (Full-length)
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9398
Gene Symbol	CD101
UniProt ID	IGSF2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---