

Anti-TOMM40 antibody [SRM14756] (STJ11109486)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IP/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:2000-1:10000 IP, 0.5 Mu g-4 Mu g antibody for 400 Mu g-600 Mu g extracts of whole cells IHC-P, 1:100-1:500 IF/ICC, 1:100-1:500 ELISA, Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on yo
Isotype	IgG
Molecular Weight	Protein Weight 103kDa Observed Mw
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10452
Gene Symbol	TOMM40
UniProt ID	TOM40_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---