

**Anti-ENPP3/CD203c antibody [SRM14760] (STJ11109490)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IF/ICC/FC/ELISA |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                         |
| Conjugation         | Unconjugated                                                                                                                                                                                       |
| Purification        | Affinity purification                                                                                                                                                                              |
| Dilution Range      | WB, 1:1000-1:2000<br>IF/ICC, 1:200-1:800<br>FC, 1:500-1:1000<br>ELISA, Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Isotype             | IgG                                                                                                                                                                                                |
| Molecular Weight    | Protein Weight 74 kDa<br>Observed Mw 74 kDa                                                                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                                                     |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">5169</a>        |
| Gene Symbol | <a href="#">ENPP3</a>       |
| UniProt ID  | <a href="#">ENPP3_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|