

Anti-PU.1/SPI1 antibody [SRM14779] (STJ11109509)
GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/IF-P/IHC-P/FC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB, 1:1000-1:6000 IF/CC, 1:50-1:200 IF-P, 1:50-1:200 IHC-P, 1:500-1:5000 FC (intra) , 1:100-1:500 ELISA, Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements
Isotype	IgG
Molecular Weight	Protein Weight 42kDa Observed Mw 45kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6688
Gene Symbol	SPI1
UniProt ID	SPI1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---